

HighFinesse  
The Standard of Accuracy

HighFinesse Tutorial

# HighFinesse Analog PID Option (1 or 2 Channels): Laser Control

## How to ...


# ... set up the HighFinesse Analog PID Option (1 or 2 channels): Laser Control

The analog output at the wavelength meter will provide up to  $\pm 4$  V. This tutorial is intended to give you a brief overview of how to configure the HighFinesse laser control settings. The tutorial does not replace reading the manual. Make sure you have read and understood it (especially section 3.5) before you start the regulation. Setting voltage bounds incorrectly might cause damage to your laser.

Here we assume that the laser is already successfully connected with a fiber to the wavelength meter. If you have any questions about that refer to the quickstart guide “HighFinesse Wavelength Meter”.


Quick Start Guide  
HighFinesse Wavelength Meter

[https://www.highfinesse.com  
/en  
/support  
/quick-start-guide.html](https://www.highfinesse.com/en/support/quick-start-guide.html)



The screenshot shows the 'Settings' menu of the software. The menu is open, and the 'Laser Control Settings ...' option is highlighted with a red box. The menu items are as follows:

- Exposure 1
- Result unit
- Pulse
- Precision
- Modes
- Display
- Laser Control Settings ...**
- Switch Settings
- Pulse Settings ...
- COM Port Settings
- Extra Settings ...



**Laser Control**

☒ Regulation signal active Port 1 2

Modify | Altering sensitivity | Error signals | Calibration  
 Reference | Regulation & Sensitivity | Bounds | Various

Signal bounds [mV]

Minimum  Maximum

☒ Adjust reference midway (0.0 mV)  
☐ Adjust reference at  mV

Behaviour on exceeding bounds

☒ Only cut at signal bounds  
☐ Output errorvalue [mV]  
 at min.  at max.


☒ Clear integral history

☐ Maximum shot-per-shot change [mV]

☒ allow towards zero  
☐ drive immediately

Move to the **frame bounds** to enter them correctly.

3




Connect the voltage output of the wavelength meter to the input at the laser  
(minimum impedance at the laser input 10 kOhm).



5

Use the **PIDSim2 Tool** to simulate good starting parameters.

Alternatively, you can set PID parameters manually in the laser control settings/frame: **“Regulation & Sensitivity”**.



6

## Start ...



PIDSim2.exe

Type: Application

Start the PIDSim2 application  
located in the path ...

Installation Path of the  
Wavelength Meter Software  
\\Tools  
\\PIDSim2.exe

... and make sure you can measure.



7

**General**

Measurement count: 200

Time / Measurement [s]: 0.020

**Laser**

1 / Amplification: 1.00 V / pm

0 V Wavelength: 632.991075

Clicking on the **small black triangle** will enable more settings.

7a

## PIDSim2 Settings

General

Measurement count

200

▲

▼

Time / Measurement [s]

0.020

▲

▼

Response count

2

▲

▼

VOut minimum [V]

0.000

▲

▼

VOut maximum [V]

4.096

▲

▼

VOut resolution [mV]

0.500

▲

▼

Synchronize f and r

☐

Laser

1 / Amplification

+

1.00

▲

▼

V / pm

0 V Wavelength

632.991075

▲

▼

Now you can **make your settings**  
for simulation.

Measurement count:

number of points in the simulation.

Time/Masurement [s]:

get this live from the wavelength meter by clicking on the red dot and confirm by clicking on the checkmark.

## Response Count:

2 for single channel, 1 for multichannel measurements.

Set the **minimum and maximum output** according to your system (voltage bounds set in 1.).

Set the **resolution** to obtain a realistic simulation of your system.





9

Click on the blue triangles to enter the bounds **9a** (range should be smaller or equal to the bounds in 1.) and number of collection points used for calculation.

Then click on the red dot **9b** to automatically **determine the sensitivity**. Once this is determined transfer the result to the frame “Regulation” and enter it as the “Sensitivity” of the laser.

**Caution:** this will vary the output voltage, so a safe choice for the bounds is important.

### Amplification settings

Collection points per cycle	900
Upper voltage border [mV]	1250
Lower voltage border [mV]	-1250

Laser	
1 / Amplification	
+ 1.00 -	V / pm
0 V Wavelength	632.991075
Perturbation [pm]	
0	
Noise [pm]	0.000
Sudden hops [pm]	0.000

**9b** Click to determine the sensitivity ...

**9a** Click to enter the bounds ...

10

Reference Voltage

Reference Voltage [mV] 1000

x

10b

Laser

1 / Amplification

+  
-

1.00

+  
-

V / pm

•  
•  
•

0 V Wavelength 632.991075

+  
-

Perturbation [pm]

0

Noise [pm] 0.000

+  
-

Sudden hops [pm] 0.000

+  
-

10c Get the current wavelength

10a Click to enter the current voltage ...

Click on the blue triangle 10a.

Enter the current voltage 10b.

Get the corresponding wavelength by clicking on the red dot 10c.

Regulation

f(t) [nm]

= 632.9910722

☒ ta [s]
 

0.040

▲▼

P

0.25

▲▼

I

1.00

▲▼

D

0.060

▲▼

Auto TPID

☐ Clear history on range exceed

Sensitivity

+

-

1.00

▲▼

V / pm

▼

Synchronize WLM



Regulation

f(t) [nm]

= 632.9910722

☒ ta [s]
 

0.040

+

-

P

0.25

+

-

I

1.00

+

-

D

0.060

+

-

Auto TPID

☐ Clear history on range exceed

Sensitivity

+

-

1.00

+

-

V / pm

Synchronize WLM

Regulation

$f(t)$  [nm<sup>2</sup>]  
= 632,9910722

☒  $t_a$  [s]  
0.040

Auto TPID

☐ Clear history on range exceed


Sensitivity  
1.00 V / pm

Synchronize WLM

Click to delete the PID parameter


**Click to determine the PID parameters.**

Press **“Auto TPID”** to determine the PID parameters. **Finally, you can close the PIDSim2.** Now the system should be ready for a test.



For this **start the Regulation.**

You can **optimize the regulation** further by **using the LongTerm application** and **minimizing possible unwanted effects** by **altering the PID parameters**.



HighFinesse  
The Standard of Accuracy



HighFinesse GmbH  
Neckarsulmer Straße 5  
72072 Tübingen, Germany



+ 49 (0) 7071 - 53 918 0  
[info@highfinesse.com](mailto:info@highfinesse.com)  
[www.highfinesse.com](http://www.highfinesse.com)



Find further information on  
products, data sheets and  
distributors on our website