

HighFinesse
The Standard of Accuracy

HighFinesse Tutorial

Extracting data from the Wavelength Meter WS6 Standalone

Important notice

Make sure the proper version of Visual C++ Redistributable for Visual Studio for your operating system is installed. You can download it here:

<https://support.microsoft.com/en-us/help/2977003/the-latest-supported-visual-c-downloads>

Then please run the vc_redist.x64.exe (vc_redist.x86.exe for 32 bit).



This tutorial shows you how to ...

... extract data from the HighFinesse Wavelength Meter WS6 Standalone using the HighFinesse NetGUI and LongTerm application.


How to ...

... use your own code to access a HighFinesse Wavelength Meter.

HighFinesse Tutorial

HighFinesse NetAccess

1



Connect the Wavelength Meter
to the Network.

2

Connect the **USB stick** with the
HighFinesse Standalone Accessories
to your computer.

If you have lost the USB stick download
the files using the link below:

[https://
www.highfinesse-downloads.com/
download/
6mcxfh5nw96w](https://www.highfinesse-downloads.com/download/6mcxfh5nw96w)



Find the RAR archive
HighFinesse
Standalone Accessories
on the **USB stick**.

The screenshot shows a file manager interface with a 'Tools' menu open. The 'Extract' option is selected, and its sub-menu is displayed. The sub-menu contains three options: 'Extract all' with the keyboard shortcut 'Ctrl+L', 'Extract displayed object(s)' with 'Ctrl+S', and 'Extract selected object(s)' with 'Ctrl+E'. The main menu also includes options like 'Add', 'Convert', 'Text encoding', 'More', 'Enter password / keyfile', 'Navigation', 'Sort by', 'Select all (Ctrl+A)', 'Select...', 'Preview with...', 'Open in a new tab', 'File manager', 'Web search', 'Explore path' (Alt+F7), and 'Properties' (Ctrl+Alt+F7).

Option	Keyboard Shortcut
Extract all	Ctrl+L
Extract displayed object(s)	Ctrl+S
Extract selected object(s)	Ctrl+E


Extract the RAR archive
HighFinesse
Standalone Accessories.



HighFinesse Standalone Accessories


Open the extracted folder
HighFinesse
Standalone Accessories.

6




Copy the folder
NetGUI for generating the correct wlmData.ini to your computer
(for WS6 Standalone only).

7



Copy the folder
Example Longterm to your computer.


8



Open the folder **NetGUI for generating the correct wlmData.ini**
(for WS6 Standalone only).

9

Open ...




Start the **NetGUI** by opening **wlm_ws6net** to test the connection.

In order to use the program UPD has to be allowed in your network.
The NetGui can be used to adjust the measurement settings and monitor the interferograms.

10

Close ...



Close the **NetGUI**.


11



Result:

The **NetGUI** has automatically overwritten the file **wlmData.ini** located in the same folder.

12



13

wlmData - Notepad


File Edit Format View Help

```
[wlm_ws6net.dat]
version=4
address=192.168.13.161
port=7171
port2=7172
offload=1
```

The file shows you the **IP address** of the standalone instrument in your network.


In this example it is **192.168.13.161**

14




Copy the **IP address** to your clipboard and close the folder and the **wlmData.ini**.

15



Open the folder **Example LongTerm ...**

16



... and open the file **wlmData.ini** with an texteditor.

17

```

wlmData - Notepad
File Edit Format View Help
; #####
; @file wlmData.ini
; @brief Example configuration file for HighFinesse NetAccess WLM/LSA network solution
; client side library
; @date: 2021.03.10
; @version: 0.1
; #####

; wlmData.ini example scenario 1 configuration file
; #####

[default]                ; Default settings (it match with all measurement applications)
version = 4                ; IPv4
address = 192.168.13.183    ; Instrument server IP address
port = 7171                ; Set/Get TCP Port number
port2 = 7172               ; CallbackProc/Ex TCP Port number
offload = 1                ; ConvertUnit / ConvertDeltaUnit functions
                           ; network offload (1=On, 0=Off)
loglevel = 3               ; Loglevel: Warning
errormode = 9              ; Error signaling: Log+EXIT 1 + 8 = 9

; wlmData.ini example scenario 2 configuration file
; #####

; [LongTerm1.exe]          ; Configuration section for LongTerm1.exe
; version = 4              ; IPv4
; address = 192.168.10.2    ; Instrument server IP address
; port = 7171              ; Set/Get TCP Port number
; port2 = 7172             ; CallbackProc/Ex TCP Port number
; offload = 1              ; ConvertUnit / ConvertDeltaUnit functions network offload (1=On, 0=Off)
;
; [LongTerm2.exe]          ; Configuration section for LongTerm2.exe

```

Find the **factory set IP address** in this file.

18



Paste ... from clipboard

Replace the **factory set IP address** by the address in your clipboard.

19




Save ... and close

Save the file **wlmData.ini**
... and close this file.

20

Open ...



Start the **LongTerm.exe**.

This program allows you to **log**
the wavelength and additional
measurement data.

See also: **Introduction LongTerm.pdf**

It also allows you to
save the data as an ASCII file.


More detailed description
of the application, see:

How to ...

... control the Wavelength Meter
with your own application via
the network.

HighFinesse Tutorial

NetAccess_UserGuide.pdf



HighFinesse
The Standard of Accuracy



HighFinesse GmbH
Neckarsulmer Straße 5
72072 Tübingen, Germany



+ 49 (0) 7071 - 53 918 0
info@highfinesse.com
www.highfinesse.com



Find further information on
products, data sheets and
distributors on our website