

HighFinesse  
The Standard of Accuracy

HighFinesse Tutorial

# Control lasers with the HighFinesse Multichannel Laser Controller (MCLC)

## How to ...

... set up the HighFinesse Multichannel Laser Controller (MCLC)

This tutorial is intended to give you a brief overview of how to configure the HighFinesse Multichannel Laser Controller. The tutorial does not replace reading the manual. Make sure you have read and understood it (especially section 3.1 and 3.2) before you start the regulation. Setting voltage bounds incorrectly might cause damage to your laser.

Here we assume that the laser is already successfully connected with a fiber to the MCLC, the MCLC is connected to the wavemeter and laser control outputs are connected to the respective laser controllers. For more information check the quickstart guide. If you have any questions about that refer to the quickstart guide “HighFinesse Multichannel Laser Controller”.


# Quick Start Guide

## HighFinesse Wavelength Meter


<https://www.highfinesse.com/en/support/quick-start-guide.html>



1



- 1a Connect the laser to the MCLC input. —
- 1b Connect the MCLC output fiber to the wavelength meter. —
- 1c Connect the MCLC laser control output to the laser controlling unit. —
- 1d Connect the wavelength meter and MCLC to the computer via the USB cable and Install the wavelength meter software. —



The screenshot shows the 'Settings' menu on the left and the 'Laser Control' window on the right. In the 'Settings' menu, the 'Settings' option is highlighted with a red box, and the 'Laser Control Settings ...' option is also highlighted with a red box. The 'Laser Control' window shows the 'Regulation & Sensitivity' tab selected, with the 'Bounds' sub-tab highlighted. The 'Reference' section shows a course value of 632.9910722 nm. The graph below shows a horizontal blue line at 632.99107 nm on a plot of intensity versus time (t[s]).



**Laser Control**

☐ Regulation signal active

Port 1 2 3 4 5 6 7 8

Modify Altering sensitivity Errorsignals Calibration  
 Reference Regulation & Sensitivity Bounds Various

Signal bounds [mV]

Minimum -10000 Maximum 10000

☒ Adjust reference midway (0.0 mV)  
☐ Adjust reference at 0 mV

Behaviour on exceeding bounds

☒ Only cut at signal bounds  
☐ Output errorvalue [mV]  
 at min. -2800 at max. -3200


☒ Clear integral history

☐ Maximum shot-per-shot change [mV]  
 10000 ☒ allow towards zero  
☐ drive immediately

Move to the **frame Bounds** to enter the minimum and maximum value (correctly).

## 4


## Using multiple channels simultaneously



Choose the port where the voltage will be put out by clicking on the **black** numbers.

4a


**Click to  
assign ...**



Then click on one of the colored numbers to assign the switch signal to the port.

**4b**

**Click again to unassign ...**




You can **unassign** it by clicking on the **same number again**. In the example the Switch signal 5 is assigned to port 2.

5

Use the **PIDSim2 Tool** to simulate good starting parameters.

Alternatively, you can set PID parameters manually in the laser control settings/frame: **“Regulation & Sensitivity”**.



6

## Start ...



PIDSim2.exe


Type: Application

Start the PIDSim2 application  
located in the path ...

Installation Path of the  
Wavelength Meter Software  
\\Tools  
\\PIDSim2.exe

... and make sure you can measure.

7



7a



95

number of points in the simulation.

get this live from the wavelength meter by clicking on the red dot and confirm by clicking on the checkmark.

2 for single channel, 1 for multichannel measurements.

according to your system (voltage bounds set in 1.).

simulation of your system.



Regulation

$f(t)$  [nm]

= 632,9910722

☒  $t_a$  [s]      P 0.25

0.040      I 1.00

Auto TPID      D 0.060

☐ Clear history on range exceed

Sensitivity

+ 1.00 V / nm

Synchronize WLM

## Connection with WLM

In which direction shall be synchronized first?

WLM --> PIDSim

PIDSim2 --> WLM

Choose ...

Choose to **synchronize** the  
PIDSim2 in the **section regulation**.

You can **alternatively also transfer all settings** you have made from the wavelength meter to the PIDsim2 tool.

Regulation

$f(t)$  [nm]

= 632,9910722

☒ ta [s] P 0.20

0.040 I 1.00

Auto TPID D 0.060

☐ Clear history on range exceed

Sensitivity 1.00

Port 1

Release WLM



9

Click on the blue triangles to enter the bounds **9a** (range should be smaller or equal to the bounds in 1.) and number of collection points used for calculation.

Then click on the red dot **9b** to automatically **determine the sensitivity**. Once this is determined transfer the result to the frame “Regulation” and enter it as the “Sensitivity” of the laser.

**Caution:** this will vary the output voltage, so a safe choice for the bounds is important.

### Amplification settings

Collection points per cycle	900
Upper voltage border [mV]	1250
Lower voltage border [mV]	-1250

Laser	
1 / Amplification	
+ 1.00 -	V / pm
0 V Wavelength	632.991075
Perturbation [pm]	
0	
Noise [pm]	0.000
Sudden hops [pm]	0.000

**9b** Click to determine the sensitivity ...

**9a** Click to enter the bounds ...

10

Reference Voltage x

Reference Voltage [mV] 1000 10b

Laser

1 / Amplification

+ 1.00 - V / pm

0 V Wavelength 632.991075 10c

Perturbation [pm]

0

Noise [pm] 0.000

Sudden hops [pm] 0.000

10c Get the current wavelength

10a Click to enter the current voltage ...

Click on the blue triangle 10a.

Enter the current voltage 10b.

Get the corresponding wavelength by clicking on the red dot 10c.

Regulation

f(t) [nm]

= 632.9910722

☒ ta [s]
 

0.040

▲

▼

P

0.25

▲

▼

I

1.00

▲

▼

D

0.060

▲

▼

Auto TPID

☐ Clear history on range exceed

Sensitivity

+

-

1.00

▲

▼

V / pm

▼

Synchronize WLM



Regulation

f(t) [nm]

= 632.9910722

☒ ta [s]
 

0.040

+

-

P

0.25

+

-

I

1.00

+

-

D

0.060

+

-

Auto TPID

☐ Clear history on range exceed

Sensitivity

+

-

1.00

+

-

V / pm

▼

Synchronize WLM

Regulation

$f(t)$  [nm<sup>2</sup>]  
= 632,9910722

☒  $t_a$  [s]  
0.040

P 0.25

I 1.00

D 0.060

Auto TPID

☐ Clear history on range exceed


Sensitivity  
1.00 V / pm

Synchronize WLM

Click to delete the PID parameter


**Click to determine the PID parameters.**

Press **“Auto TPID”** to determine the PID parameters. **Finally, you can close the PIDSim2.** Now the system should be ready for a test.



For this **start the Regulation.**

You can **optimize the regulation** further by using the LongTerm application and minimizing possible unwanted effects by altering the PID parameters.



HighFinesse  
The Standard of Accuracy



HighFinesse GmbH  
Neckarsulmer Straße 5  
72072 Tübingen, Germany



+ 49 (0) 7071 - 53 918 0  
[info@highfinesse.com](mailto:info@highfinesse.com)  
[www.highfinesse.com](http://www.highfinesse.com)



Find further information on  
products, data sheets and  
distributors on our website