



Combination  
Calibration Source  
Wavelength Meter



Wavelength	Accuracy	Wavelength Meter	Calibration Source	Frequency stability
633 nm	30 MHz	WS7-30	SL04	± 5 MHz
	60 MHz	WS7-60	SL04	± 5 MHz
			SL02	± 25 MHz

Wavelength	Accuracy	Wavelength Meter	Calibration Source	Frequency stability
780 nm	2 MHz	WS8-2	SLR 780	± 0.5 MHz
	10 MHz	WS8-10	SLR 780	± 0.5 MHz
	30 MHz	WS7-30	SLR 780	± 0.5 MHz
	60 MHz	WS7-60	SLR 780	± 0.5 MHz

Wavelength	Accuracy	Wavelength Meter	Calibration Source	Frequency stability
852 nm	2 MHz	WS8-2	SLR 852	± 2 MHz
	10 MHz	WS8-10	SLR 852	± 2 MHz
	30 MHz	WS7-30	SLR 852	± 2 MHz
	60 MHz	WS7-60	SLR 852	± 2 MHz

Wavelength	Accuracy	Wavelength Meter	Calibration Source	Frequency stability
1532 nm	2 MHz	WS8-2	LFR-1532	± 2 MHz
	10 MHz	WS8-10	LFR-1532	± 2 MHz
	30 MHz	WS7-30	LFR-1532	± 2 MHz
	60 MHz	WS7-60	LFR-1532	± 2 MHz



HighFinesse GmbH  
Wöhrdstraße 4  
72072 Tübingen, Germany



T + 49 (0) 7071 - 539180  
F + 49 (0) 7071 - 5391899  
M info@highfinesse.com



Additional information  
and distributors:  
www.highfinesse.com