

Available Measurement Ranges

PCF	Singlemode, Range 370 – 2000 nm
-----	---------------------------------

Switch


Available channels	4, 8
Guaranteed transmission	15 % @ 633 nm
Return loss	< 40 dB
Crosstalk	< –55 dB
Connectors	FC/APC
Switching time	typ. 12 ms

Control Outputs

Available channels	4, 8
Resolution	16 bit
Voltage span	highZ –10 V .. +10 V
Connector	Lemo 00 and BNC adapters included in the scope of delivery

General

Supply voltage	9 – 24 V
Power consumption	max. 4 W
Operating temperature	0 – 40° C, non condensing
Storage temperature	–40° – 80° C
Housing dimensions	170 mm × 175 mm × 55 mm (without fiber)
Standard interface	USB 2.0



Multi Channel Laser Controller

MCLC



HighFinesse

The Standard of Accuracy

Further Information

For further technical information, application examples, diagrams
and for customisation of singlemode switches please contact:

HighFinesse Service

service@highfinesse.de



HighFinesse GmbH
Neckarsulmer Straße 5
72072 Tübingen, Germany



T +49 (0) 7071 - 53 918 0
F +49 (0) 7071 - 53 918 99
M info@highfinesse.com



Additional information
and distributors:
www.highfinesse.com