



Switches

Photonic Crystal Switch Rack




HighFinesse
The Standard of Accuracy

Measurement Range

370 – 2000 nm

Technical Specs

Available channels	4, 8
Guaranteed transmission	15 % @ 633 nm
Return loss	< 40 dB
Crosstalk	< -55 dB
Switching time	typ. 12 ms
Supply voltage	powered by wavelength meter (WR Series)
Operating temperature	0 – 60° C
Storage temperature	-40° – 80° C
Operating Relative Humidity	35 – 55 %
Dimensions	360 mm × 483 mm × 43.5 mm
Weight	4.0 kg
Standard interface	TTL, USB for service only





Switches

Photonic Crystal Switch Rack



HighFinesse
The Standard of Accuracy

Further Information

For further technical information, application examples, diagrams and for customisation of photonic crystal switch racks please contact:

HighFinesse Service

service@highfinesse.de



HighFinesse GmbH
Neckarsulmer Straße 5
72072 Tübingen, Germany



T +49 (0) 7071 - 53 918 0
F +49 (0) 7071 - 53 918 99
M info@highfinesse.com



Additional information
and distributors:
www.highfinesse.com